

Invacare[®] Nutron[™]R50∟x Power Wheelchair

Invacare Nutron R50LX Power Wheelchair

s it possible for one power chair to be both stylish and powerful? The Nutron[™] R50L× is both. Its compact design and dependable rear-wheel drive make it easy to maneuver around obstacles inside and outside the home. And when the consumer wants to travel, the R50L× disassembles and folds for easy transport.



Specifications

Overall dimensions	
Length*	40.75"
Widths (excluding joystick)	20" - 26"
Heights	
17.75" seat-to-floor	35" - 38"
19.75" seat-to-floor	37" - 40"
Seat dimensions	
Widths	14" - 20"
Depths	16" - 18"
Back heights	16" - 19"
Adjustable-angle back	-5° to +20° (5° increments)
Seat angle	3° fixed
Weight capacity	250 lb.
Wheel and caster sizes	
Front casters (standard)	8" × 2"
Front casters (optional)	6" x 2"
Drive wheels	12.5" x 2.25"
Maximum speed	4 mph
Incline capability	9 °
Turning radius**	
Front (footrest in place)	33"
Rear (anti-tippers in place)	12"
Weight	
With 22NF batteries	I 55 lb.
Without batteries and front riggings	85 lb.
Ground clearance	3.5"
Warranties	
Base frame	60 months
Seat frame	36 months
Electronics	12 months
Motor/gearbox	12 months
* M I D M (02.6	

* Measured with PW93 front riggings and anti-tippers with 16" seat depth.

** Measured with 16"W x 16"D seat.



The R50LX can be disassembled and folded for easy transport.

Features

- Rear-wheel drive and compact size provide optimum maneuverability and accessibility
- Oversize, high-strength tubular frame with bold styling
- Seat-to-floor heights as low as 17.75" permit easier transfers
- Invacare MKIV[™]-RII LP electronics for precise driver control
- Reinforced crossbrace design folds to 13" for transport
- Back angle and height adjust for maximum comfort and support
- Telescoping front riggings allow optimal foot placement
- Flat-free wheel inserts mean worry-free operation
- Powerful 350-watt, high-torque motors ensure dependable performance
- Regenerative braking and automatic motor locks for smooth stopping
- Nylon upholstery adds style and comfort
- 8-amp charger provides maximum power
- Waterfall arm pads offer user support
- Optional fold-down back permits maximum portability
- 250-lb. weight capacity
- Chair is available in 5 frame colors



Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada 570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare is a registered trademark and Yes, you can., Nutron and MKIV are trademarks of Invacare Corporation. Specifications are subject to change <u>without</u> notification.